

Package: LEEF.measurement.flowcytometer (via r-universe)

September 16, 2024

Type Package

Title What the Package Does (Title Case)

Version 0.8.13

Maintainer Rainer M. Krug <Rainer.Krug@uzh.ch>

BugReports <https://github.com/LEEF-UZH/LEEF.measurement.flowcytometer/issues>

URL <https://github.com/LEEF-UZH/LEEF.measurement.flowcytometer>

Description More about what it does (maybe more than one line) Use
four spaces when indenting paragraphs within the Description.

License MIT + file LICENSE

Depends R (>= 3.5.0)

Imports R.utils, flowCore (>= 2.0), tiff, yaml, plyr, tidyr, loggit,
magrittr

Encoding UTF-8

LazyData true

RoxygenNote 7.1.2

Suggests testthat, openssl

Repository <https://leef-uzh.r-universe.dev>

RemoteUrl <https://github.com/LEEF-UZH/LEEF.measurement.flowcytometer>

RemoteRef master

RemoteSha bfa9bb3121db3e16150c657dfc91a311f6963b41

Contents

add_new_data	2
extractor_flowcytometer	2
extractor_flowcytometer_gating	3
extractor_flowcytometer_preparation	3
load_parameter	4

par_template	4
pre_processor_flowcytometer	5
pre_processor_flowcytometer_c6	5
print_parameter	6
raw_data_ok	6
register	7
save_parameter	7

Index	8
--------------	----------

add_new_data	<i>Check if data in input folder is OK and move to raw data folder</i>
--------------	--

Description

Check if data in input folder is OK and move to raw data folder

Usage

```
add_new_data(input, output)
```

Arguments

input	The folder, where a folder flowcytometer is located which contains the new files.
output	A folder, which contains a subfolder called flowcytometer, i.e. the usually the raw data folder, into which the files will be moved to.

Value

a list which contains the individual results for each file. TRUE if moved, FALSE if an error occurred. Details of the error re in the error files in the input/flowcytometer directory.

extractor_flowcytometer	<i>Extractor flowcytometer data</i>
-------------------------	-------------------------------------

Description

This function is extracting data to be added to the database (and therefore make accessible for further analysis and forecasting) from .fcs files.

Usage

```
extractor_flowcytometer(input, output)
```

Arguments

input	directory from which to read the data
output	directory to which to write the data

Value

invisibly TRUE when completed successful

`extractor_flowcytometer_gating`
Extractor flowcytometer data

Description

This function is extracting data to be added to the database (and therefore make accessible for further analysis and forecasting) from .fcs files.

Usage

```
extractor_flowcytometer_gating(input, output)
```

Arguments

input	directory from which to read the data
output	directory to which to write the data

Value

invisibly TRUE when completed successful

`extractor_flowcytometer_preparation`
Extractor flowcytometer preparation data

Description

This function is extracting data to be added to the database (and therefore make accessible for further analysis and forecasting) from .fcs files.

Usage

```
extractor_flowcytometer_preparation(input, output)
```

Arguments

input	directory from which to read the data
output	directory to which to write the data

Value

invisibly TRUE when completed successful

load_parameter	<i>Load parameter from file</i>
----------------	---------------------------------

Description

Load parameter from file

Usage

```
load_parameter(file = "parameter.yaml")
```

Arguments

file	name of parameter file
------	------------------------

Value

invisibly TRUE

par_template	<i>Template function to assign value to parameter in the package wide cache</i>
--------------	---

Description

assign the function to a new cvariable and the name of the function will be used for the parameter name. e.g:

- `fps <- par_template`

Usage

```
par_template(value)
```

Arguments

value	if missing, the value of the parameter will be returned, NULL if the parameter does not exist; if specified, the parameter will be set to the value
-------	---

Value

the (new) value of the argument

pre_processor_flowcytometer
Preprocessor flowcytometer data

Description

Just copy all files

Usage

pre_processor_flowcytometer(input, output)

Arguments

input	directory from which to read the data
output	directory to which to write the data

Value

invisibly TRUE when completed successful

pre_processor_flowcytometer_c6
Preprocessor flowcytometer data FOR c6 FILES!!!

Description

Convert all .c6 files in flowcytometrie folder to .fcs files in output folder

Usage

pre_processor_flowcytometer_c6(input, output)

Arguments

input	directory from which to read the data
output	directory to which to write the data

Value

invisibly TRUE when completed successful

print_parameter	<i>Print the flowcytometer parameter</i>
-----------------	--

Description

Print the flowcytometer parameter

Usage

```
print_parameter(print_as_yaml = TRUE, echo = TRUE)
```

Arguments

print_as_yaml	Print in yaml formatted text; ~ stands for NULL
echo	if TRUE print the parameter, if FALSE just return them as list

Value

invisible returns list of parameter for further processing

raw_data_ok	<i>Check if data in raw data folder is OK</i>
-------------	---

Description

Check if data in raw data folder is OK

Usage

```
raw_data_ok(input)
```

Arguments

input	raw data folder containing flowcytometer data, i.e usually is some/path/flowcytometer
-------	---

Value

TRUE if ok, FALSE or list of problems if not

Examples

```
## Not run:
raw_data_ok()

## End(Not run)
```

register	<i>Register the processing of flowcytometer data in the LEEF package</i>
----------	--

Description

Register the processing of flowcytometer data in the LEEF package

Usage

```
register()
```

Value

invisibly TRUE when completed successful

save_parameter	<i>Save parameter into .yaml file</i>
----------------	---------------------------------------

Description

Save parameter into .yaml file

Usage

```
save_parameter(file = "parameter.yaml")
```

Arguments

file	name of parameter file
------	------------------------

Value

invisibly TRUE

Index

`add_new_data`, [2](#)

`extractor_flowcytometer`, [2](#)

`extractor_flowcytometer_gating`, [3](#)

`extractor_flowcytometer_preparation`, [3](#)

`load_parameter`, [4](#)

`par_template`, [4](#)

`pre_processor_flowcytometer`, [5](#)

`pre_processor_flowcytometer_c6`, [5](#)

`print_parameter`, [6](#)

`raw_data_ok`, [6](#)

`register`, [7](#)

`save_parameter`, [7](#)